

HEMOSTASIS

Five level, liquid-stable products specifically designed to challenge the Siemens Sysmex system's full reportable range and eliminate manual dilutions.



VALIDATE® D-DIMER

ORDER NO. 902se **CE**

Typical Recovery:
0.3 - 3.86 mg/L FEU

Matrix:
Human Plasma

Form:
3.0 mL, Liquid-Frozen

Storage:
-10 to -25 C

Open-Vial:
4 freeze-thaws through expiration

VALIDATE® HEPARIN

ORDER NO. 903se **CE**

Typical Recovery:
0.13 - 1.29 IU/mL

Matrix:
Human Plasma

Form:
3.0 mL, Liquid-Frozen

Storage:
-10 to -25 C

Open-Vial:
2 freeze-thaws through expiration

VALIDATE® FIBRINOGEN

ORDER NO. 904se **CE**

Typical Recovery:
94 - 432 mg/dL

Matrix:
Human Plasma

Form:
2.0 mL, Liquid-Frozen

Storage:
-10 to -25 C

Open-Vial:
Single-use, 6 hours

EASY • FAST • EFFICIENT



Liquid, Ready-to-Use

Increase productivity and avoid the time, error and inconvenience of reconstitution. Supplied in dropper bottles, simply bring to room temperature and dispense into the sample cup, reducing the risk of contamination.



Open-Vial Stability

Freeze-thaw cycles allow for multiple uses - cap the bottle, return to storage and use again. For Heparin and Fibrinogen, order once and receive a second, no charge shipment, in 6 months.



Troubleshooting

VALIDATE is ideal for use after preventative maintenance or to troubleshoot across the full reportable range.



Linearity & Peer Groups

5 Level products with free data reduction services offers instant linearity results via MSDRx with ready access to Peer Statistics.



Quality

Products are U.S. FDA cleared or listed and are CE marked. Our facility is registered with the U.S. FDA and our QMS is certified to ISO 13485.



Accreditation

For use at installation, or during the laboratory's routine 6-month testing interval, rely on VALIDATE to satisfy ISO 15189, CAP, JCI and local regulatory requirements, including AMR validation.